

IT PORTA STOCK
EN PORTA STOCK
DE STOCK-HALTER
FR PORTA STOCK
ES PORTA STOCK
PT PORTA STOCK

IT PORTA STOCK CON RUOTE
EN PORTA STOCK WITH WHEELS
DE STOCK-HALTER MIT RÄDERN
FR PORTA STOCK AVEC ROUES
ES PORTA STOCK CON RUEDAS
PT PORTA STOCK COM RODAS

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IT

Instructions for use PORTA STOCK – SUPPORTING FRAME FOR FIELD STRETCHERS

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EN

1. MODELS	7
2. INTENDED USE	7
2.1 INTENDED USE	7
2.2 TARGET PATIENTS	7
2.3 PATIENT SELECTION CRITERIA	7
2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS	7
2.5 USERS AND INSTALLERS	7
3. REFERENCE STANDARDS	7
4. INTRODUCTION	7
4.1 DEVICE LABELLING AND TRACEABILITY	7
4.2 SYMBOLS	7
5. WARNINGS/DANGERS	7
6. SPECIFIC WARNINGS	8
7. RESIDUAL RISK	8
8. TECHNICAL DATA AND COMPONENTS	8
9. PROPER USE	8
10. CLEANING AND MAINTENANCE	9
10.1 CLEANING	9
10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL	9
10.3 LIFE SPAN	9
11. TROUBLESHOOTING TABLE	9
12. ACCESSORIES	9
13. SPARE PARTS	9
14. DISPOSAL	9

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1. MODELS

The following basic models may be subject to implementation or change without notice.

- PORTA STOCK
- PORTA STOCK WITH WHEELS

2. INTENDED USE

2.1 INTENDED USE

Portastock is a support for field stretchers that allows the stretcher to be kept raised from the ground, facilitating small movements and access to the patient by the doctor. Available in versions with and without wheels, it is compatible with the Stock series stretchers and is also suitable for Spencer stretchers with a maximum width of 575 mm and a minimum length of 1250 mm (not including handles). Unsupervised parking is not permitted.

2.2 TARGET PATIENTS

The target patients are those for whom use of the field stretcher is foreseen.

2.3 PATIENT SELECTION CRITERIA

The selection criteria are those foreseen for use of the field stretcher.

2.4 CONTRAINDICATIONS AND UNWANTED SIDE EFFECTS

No particular contraindications or side effects are known with relation to use of the device, as long as it is used in accordance with the user manual.

2.5 USERS AND INSTALLERS

The intended users are those foreseen for use of field stretchers.

These devices are not intended for lay people.

Operators must be able to provide the necessary patient care.

3. REFERENCE STANDARDS

REFERENCE	DOCUMENT TITLE
EU Regulation 2017/745	EU Regulation on Medical Devices

4. INTRODUCTION

4.1 DEVICE LABELLING AND TRACEABILITY

Each device is provided with a label, placed on the device itself and/or on the packaging, which contains the Manufacturer's identification data, product, CE marking, serial number (SN) or lot number (LOT). This must never be removed or covered.

EU Regulation 2017/745 requires manufacturers and distributors of medical devices to keep track of their location.

4.2 SYMBOLS

Symbol	Meaning	Symbol	Meaning
	Device in compliance with EU Regulation 2017/745		See the user manual.
	Medical device		Lot number
	Manufacturer		Product code
	Date of manufacture		Danger – Indicates a hazardous situation that may result in a situation directly related to serious injury or death.
	Unique Device Identifier		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner (only for USA Market)

(01)08057711230006(11)200626(21)1234567890

Production identification
Alphanumeric code that identifies the production units of the device, composed of:
(01)0805771123 company prefix
000 progressive GS1
6 control number
(11)200626 date of production (YYMMDD)
(21) 1234567890 lot/SN

4.3 WARRANTY AND SERVICE

Spencer Italia S.r.l. guarantees that products are free from defects for a period of **one year from the date of purchase**.

Spencer Customer Service tel. +39 0521 541154, fax +39 0521 541222, e-mail service@spencer.it.

Warranty and service conditions are available at <http://support.spencer.it>.

5. WARNINGS/DANGERS

Product features

Use of the product for any purpose other than that described in the User Manual is prohibited.

- The product must not be tampered with or modified without the manufacturer's authorisation.
- Avoid contact with sharp or abrasive objects.
- Operating temperature: from -5°C to + 50°C.
- Storage temperature: -10°C to +60°C.

General warnings for medical devices

- Do not use if the device or parts of it are punctured, torn, frayed, or excessively worn.
- Do not alter or modify the device arbitrarily, as doing so could result in unpredictable operation and damage to the patient or rescuers and shall void the manufacturer's warranty and release the manufacturer from all liability.
- Participate in safety checks on products placed on the market, transmitting information regarding product risks to the Manufacturer as well as to the Competent Authorities for their respective actions.

IT

EN

DE

FR

ES

PT

6. SPECIFIC WARNINGS

To use the product, you must also have read, understood and carefully follow all the instructions in the user manual.

- ⚠️ **Never** apply a weight greater than the maximum capacity of the device. When determining the applicable weight, the operator should consider the weight of the patient, device, equipment, accessories and any other object attached to the stretcher.
- ⚠️ **Perform simulations of use with a stretcher and a patient simulating load and accessories before putting the device into service.**
- Never leave the patient on the device unsupervised, as they could get injured.
- Do not drag the device.
- Hold the stretcher firmly on the Porta Stock if the patient is lying on it.

7. RESIDUAL RISK

No residual risks, or rather risks that could arise despite compliance with all warnings in this user manual, have been identified.


8. TECHNICAL DATA AND COMPONENTS

Component ID	Description
1	Leg Made of steel, it is the load-bearing component that provides support on the ground.
2	Support crossbars Made of steel, these are the parts on which the stretcher must be placed.
3	Tubular joint Made of steel, it holds the structure together and allows it to be folded.
4	Leg locking compass Made of steel, it allows the device to be locked in the open position.
5	Safety belt Made of nylon, it must be used to keep the locking compass locked, ensuring that the device remains in the open position. Only present on the version with wheels.
6	Wheels If present, they are made of rubber and nylon, allowing the device with the stretcher on top to be moved on flat surfaces. Two opposing wheels are equipped with a parking brake.


	PORTA STOCK	PORTA STOCK WITH WHEELS
Length when open	1325 ± 10 mm	1325 ± 10 mm
Width	640 ± 10 mm	640 ± 10 mm
Height when open	750 ± 10 mm	750 ± 10 mm
Overall dimensions when folded	870x640x130 ± 10 mm	975x640x130 ± 10 mm
Wheels	/	Ø 150mm
Weight	10 ± 0,5 kg	11 ± 0,5 kg
Capacity	150 kg	150 kg
Frame material	Steel	Steel

9. PROPER USE


Find a level area where you can place the device to ensure its stability.
To open the device:

Spread the two legs apart.


Press on the tubular joint until the two halves are collinear, then press the locking compasses until they lock.


For the model fitted with wheels, the belt on the upper part of the frame must be tightened as much as possible.


To close the device, push the locking compasses in the opposite direction to the locking direction, lift the tubular joint and bring the legs together until the fully closed position is reached.

After making sure that the support has been correctly opened and locked, the stretcher can be placed on the support, taking care to position it centred on the frame, symmetrically with respect to the support points.

Make sure that placement is stable.

To avoid accidental falls, **DO NOT ALLOW THE PATIENT TO GET ON OR OFF the stretcher when it is on the frame.**

We recommend using the parking brakes in order to provide more stability and facilitate the transfer on the wheeled version.

Simply press down the red lever on the bracket to engage them.

The brake is on two opposing wheels to provide greater stability.

When the brake is no longer needed, release it by returning the red lever to its original position.

10. CLEANING AND MAINTENANCE

10.1 CLEANING

Failure to carry out the correct cleaning operations could increase the risk of cross-infection due to presence of body fluids and/or residues.

The operator must wear suitable personal protective equipment, such as gloves, goggles, etc. during all checking and cleaning operations.

Clean the exposed parts with water and delicate soap. **Never use solvents or stain removers.**

Do not use cleaning agents containing chlorine / sodium hypochlorite as corrosion of components may occur.

Rinse thoroughly with lukewarm water, making sure you have removed all traces of soap, which may deteriorate or compromise conditions and durability. **Avoid using high pressure water**, as it penetrates the joints and removes lubricant, creating the risk of corrosion on components. Let dry completely before storing. Drying after washing or after use in a wet environment must be natural and not forced. Do not use flames or other direct heat sources.

If disinfecting, use products that do not have a solvent or corrosive action on materials constituting the device. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators.

10.2 ROUTINE MAINTENANCE AND PERIODIC OVERHAUL

Routine maintenance is the responsibility of the user, who must carry out checks before and after each use and in any case at least every 3 months:

- General functionality of the device
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections)
- Fulfilment of the requirements of the user manual in section Warnings and Specific Warnings.
- Fulfilment of the requirements of the manual in section on Proper use.

No periodic overhaul is foreseen for the device.

10.3 LIFE SPAN

The device, if used as described in the following instructions, has a life span of 10 years from the date of purchase.

Spencer Italia S.r.l. will accept no responsibility for incorrect operation or damage caused by the use of devices that have not been overhauled by the Manufacturer or authorized centre, or that have exceeded the maximum allowable life span.

11. TROUBLESHOOTING TABLE

PROBLEM	CAUSE	REMEDY
Pushing the compasses does not lock the legs.	The compasses have been forced and the stop has been broken.	Remove the device from service and contact customer service.
Damage to the structure	Improper use or exceeded the maximum capacity	Remove the device from service and contact customer service.
It is not possible to hook the buckles of the safety belt.	Possible incorrect buckle orientation or buckle failure	Check the correct orientation of the buckle and its conditions. If damaged, immediately remove the device from service and contact the manufacturer.

If a problem or fault is detected that does not correspond to the above, please contact Spencer Italia srl customer care service.

12. ACCESSORIES

There are no accessories for these devices.

13. SPARE PARTS

There are no spare parts available for these devices.

14. DISPOSAL

When devices and their accessories are no longer suitable for use, they can be disposed of as normal municipal solid waste if they have not been contaminated by special agents. Otherwise, follow the regulations in force regarding disposal.

Warning

The information contained in this document is subject to change without notice and is to be intended as a commitment by Spencer Italia S.r.l. subject to change. The images are included as examples and may vary slightly from the actual device.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of the document may be photocopied, reproduced or translated into another language without prior written consent from Spencer Italia S.r.l.

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